

NMR New User Workshop

Second session

Rui Huang

Notes reproduced from Dr. Weiguo Hu

Link:

ruihuang@umass.edu

weiguoh@polysci.umass.edu


Booklet:

<https://udrive.oit.umass.edu/weiguoh/nmrhandouts/Introduction-NMR.pdf>

Questions and Hands-On Practices

- What are some safety concerns in the NMR lab?
- What are some best practices when preparing NMR tubes for experiments?
- What are the purposes of atma, lock, shim, and rga?
- What are some necessary spectral processing steps to take before integrating the peaks?

^{13}C NMR



- edc and select CARBON1 in experiment (attach the screenshot)
- *atma*
- Sensitivity of ^{13}C : ^1H
 - The more sample amount, the better the s/n will be
 - Using line broadening (lb) to enhance s/n
- Every carbon has one single peak
 - Because a “decoupling” technique is used to remove influence from ^1H
- Integration is less reliable than ^1H

^{13}C integration can be very useful if used wisely


<https://blogs.umass.edu/weiguoh/>



Decoupling

In ^{13}C experiment, ^{13}C - ^1H J-coupling is usually removed by a “decoupling” pulse sequence

- Pulsing at ^1H during detection of ^{13}C signal
 - i.e. Set pulsing frequency to ^1H and detector frequency to ^{13}C
- Decoupling resolves overcrowding due to large ^1J splitting
 - The large ^1J couplings usually don't provide useful information
- Decoupling reduces ^{13}C multiplets to singlets, which improves S/N

Coupling ^{13}C ^1H







DEPT (Distortionless Enhancement by Polarization Transfer)

Discern methyl, methylene, methine, and quaternary carbons

- A ^{13}C technique, with several variations. Most popular is DEPT-135
- CH and CH₃ peaks are positive; CH₂ peaks are negative
- Carbons with no directly bonded protons do not have signals
- *Why do solvent peaks vanish in DEPT?*



2D NMR : implementation


- Conventional NMR spectra (one-dimensional spectra) are plots of intensity vs. frequency
- In two-dimensional spectroscopy intensity is plotted as a function of two frequencies.

COSY : ^1H and ^1H correlation : *Connectivity of protons*

HMQC: ^1H and ^{13}C correlation : *Direct C – H Connectivity*


- **edc** to create the directory
 - Running a 2D experiment is almost as easy as running a 1D
- All the pulse sequence parameters have been “packaged” for you

COSY(Correlation Spectroscopy)



- Detects connectivity between protons that have J coupling
- • Both dimensions are ^1H
- • Vertical is always 1st dim.
- • 3 diagonal peaks
 - provide no new information
- • Cross peaks at (1.2ppm, 2.6ppm) and (2.6ppm, 1.2ppm) indicate strong J-coupling
 - Usually only 2 J and 3 J are strong enough to give cross peaks
- • You can optionally plot 1D spectrum on the side to compare with the 2D spectrum. 1D spectrum usually has higher resolution and less artifacts.
- • **Fourier transform command:** *xfb* instead of *efp*
 - *xfb* be can done at any point during the analysis
 - when good resolution is achieved the experiment can be stopped
 - using *halt* command

HMQC (Heteronuclear Multi-Quantum Correlation)



A ^{13}C - ^1H 2D spectrum

- No diagonal peaks
 - Diagonal peaks could only exist when both dimensions are the same type of nuclei
- Shows direct C-H connectivity
- Non-protonated carbons don't have a peak on the spectrum
- Beware of artifacts
 - Along the "ridges" of big peaks
 - In the middle lines of spectrum window
- Must do *atma*

NMR Facility Policies

1.Reservation

- please follow the reservation rules on top of the Calcium calendars.
- 3 violations will face suspension of their accounts for 3 weeks.

2.Internet access is limited to umass domain

3.Data storage: The NMR computer gets full quickly. You need to transfer out data and delete them from the NMR computer asap.

- pse server automatically deletes data that are > 7 days old
- Don't use flash drive; it could contaminate
- If you have too many data files, Tospin will load very slowly

- You have two folders to store your data:
 - /home/[username]. This includes your desktop. Intended for temporary storage only.
 - /opt/topspin3.2/data/[username]/nmr. You can keep stuff here for slightly longer.

NMR Facility Policies

5. Reporting of problems and incidents

- Report instrument problems to Czar if your group has one, otherwise report to me.
- Never attempt to fix problem or reboot computer
- Report incident caused by you (standard sample breakage etc.) to Dr.Hu

6. You are responsible for the act of anybody using your account (junior grad students, summer students et al)

- When you feel your trainees can comfortably run simple NMR experiments themselves, ask them to contact me to obtain their own accounts

7. Please keep the workplace clean